



SETTING THE STANDARD FOR TABLE POTATOES

Produces high yields of uniform yellow fleshed tubers with shallow eyes and superb skin finish.

Great taste and texture, a true all rounder for table use.

Dependable long-term storability.

Good overall disease resistance.

ADVICE

Do you need specific advice about seed rates, fertilisers or agronomy?
Ask your local supplier or Meijer contact.

Sensitive to Sencorex (Metribuzin) used at the full dose when applied pre-emergence.

Cure tubers thoroughly with ventilation post-harvest.

STORAGE

2°C

9 Months

CHARACTERISTICS

Maturity	6
Skin Colour	Light Yellow
Fresh Colour*	7
Tuber Shape**	RO
Eye Depth	7
Tuber Uniformity	7
Tuber Size	7.5
Freedom from Outgrades	8
Dry Matter %	21
Suitability for Cooking †	8AB
Dormancy	8

COMMENTS

1 = Negative 9 = Positive
 * 8 = Yellow 4 = White
 ** R = Round O = Oval L = Long
 † A = Firm D = Mealy

RESISTANCE

PLRV	8	Resistance to Nematodes Ro/Pa(x) = Globodera rostochien- sis/G.pallida. Pathotype(x) Ro1(=A), Ro2/Ro3(=BC), Pa2(=D), Pa3(=E) 9 = Highly Resistant 1 = Very Susceptible
A-virus	8	
X-virus	8.5	
Y-virus	7.5	
Yntn-virus	R	
Foliage Blight	4	Resistance to Wart Disease: wz 1(D1), wz 2/6 (G1/O1, wz 18 (T1) 10 = Highly Resistant 3 = Very Susceptible
Tuber Blight	7.5	
Scabies	7	
Spraing	7	
Bruising	8	
Secondary Growth	8	
Ro1	9	
Ro2 / 3	4	
Pa2	5	
Pa3	1	
wz 1 (D1)	10	
wz 2/6 (G1 / O1)	9	
wz 18 (T1)	-	